

Product Datasheet

Rat Junctional Adhesion Molecule C ELISA Kit (orb1211444)

Description	Rat Junctional Adhesion Molecule C, JAM3 ELISA Kit
Reactivity	Rat
Range	2-700ng/L
Tested Applications	ELISA
Target	JAM3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,783
Uniprot ID	Q68FQ2
Solubility (25°C)	800ng/L
Sensitivity	1.12ng/L
Entrez	315509
Expiration Date	Please enquire.