

Product Datasheet

Rat Glutaminase Kidney Isoform ELISA Kit (orb1211469)

Description	Rat Glutaminase Kidney Isoform, Mitochondria, GLS ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	GLS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	74,024
Uniprot ID	P13264
Solubility (25°C)	3200ng/L
Sensitivity	14.34ng/L
Entrez	24398
Expiration Date	Please enquire.