

Product Datasheet

Rat Farnesyl Pyrophosphate Synthase ELISA Kit (orb1211485)

Description	Rat Farnesyl Pyrophosphate Synthase, FDPS ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	FDPS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,830
Uniprot ID	P05369
Solubility (25°C)	3200ng/L
Sensitivity	7.16ng/L
Entrez	83791
Expiration Date	Please enquire.