

Product Datasheet

Rat Disks Large Homolog 4 ELISA Kit (orb1211490)

Description	Rat Disks Large Homolog 4, DLG4 ELISA Kit
Reactivity	Rat
Range	20-4500ng/L
Tested Applications	ELISA
Target	DLG4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	80,465
Uniprot ID	P31016
Solubility (25°C)	4800ng/L
Sensitivity	9.88ng/L
Entrez	29495
Expiration Date	Please enquire.