

Product Datasheet

Rat C-C Chemokine Receptor Type 2 ELISA Kit (orb1211504)

Description	Rat C-C Chemokine Receptor Type 2, CCR2 ELISA Kit
Reactivity	Rat
Range	6-1800ng/L
Tested Applications	ELISA
Target	CCR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,763
Uniprot ID	055193
Solubility (25°C)	1920ng/L
Sensitivity	3.35ng/L
Entrez	60463
Expiration Date	Please enquire.