

Product Datasheet

Rat Bcl-2-like Protein 1 ELISA Kit (orb1211511)

Description	Rat Bcl-2-like Protein 1, BCL2L1 ELISA Kit
Reactivity	Rat
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	BCL2L1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,158
Uniprot ID	P53563
Solubility (25°C)	12.8ng/ml
Sensitivity	0.012ng/ml
Entrez	24888
Expiration Date	Please enquire.