

Product Datasheet

Rat Arginase-2 ELISA Kit (orb1211514)

Description	Rat Arginase-2, ARG2 ELISA Kit
Reactivity	Rat
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	ARG2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,640
Uniprot ID	008701
Solubility (25°C)	96ng/mL
Sensitivity	0.152ng/mL
Entrez	29215
Expiration Date	Please enquire.