

Product Datasheet

Rat Apolipoprotein A-4 ELISA Kit (orb1211516)

Description	Rat Apolipoprotein A-4, APOA4 ELISA Kit
Reactivity	Rat
Range	5-1000ng/mL
Tested Applications	ELISA
Target	APOA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,456
Uniprot ID	P02651
Solubility (25°C)	1280ng/mL
Sensitivity	2.22ng/mL
Entrez	25080
Expiration Date	Please enquire.