

Product Datasheet

Rat Aminopeptidase N ELISA Kit (orb1211520)

Description	Rat Aminopeptidase N, ANPEP ELISA Kit
Reactivity	Rat
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	ANPEP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	109,449
Uniprot ID	P15684
Solubility (25°C)	48ng/mL
Sensitivity	0.523ng/mL
Entrez	81641
Expiration Date	Please enquire.