

Product Datasheet

Rat 5-Hydroxytryptamine Receptor 2A ELISA Kit (orb1211533)

Description	Rat 5-Hydroxytryptamine Receptor 2A, HTR2A ELISA Kit
Reactivity	Rat
Range	10-3000ng/L
Tested Applications	ELISA
Target	HTR2A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,850
Uniprot ID	P14842
Solubility (25°C)	3200ng/L
Sensitivity	6.89ng/L
Entrez	29595
Expiration Date	Please enquire.