

Product Datasheet

Rat Kininogen-1 ELISA Kit (orb1211560)

Description	Rat Kininogen-1, KNG1 ELISA Kit
Reactivity	Rat
Range	5-1500ng/ml
Tested Applications	ELISA
Target	KNG1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70933
Uniprot ID	P08934
Solubility (25°C)	1600ng/ml
Sensitivity	3.06ng/ml
Entrez	25087
Expiration Date	Please enquire.