

Product Datasheet

Rat Insulin-like Growth Factors 1 Receptor ELISA Kit (orb1211563)

Description	Rat Insulin-like Growth Factors 1 Receptor, IGF-1R ELISA Kit
Reactivity	Rat
Range	1-400ng/L
Tested Applications	ELISA
Target	IGF1R
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	155,396
Uniprot ID	P24062
Solubility (25°C)	480ng/L
Sensitivity	0.44ng/L
Entrez	25718
Expiration Date	Please enquire.