

Product Datasheet

Rat Lactoylglutathione Lyase ELISA Kit (orb1211580)

Description	Rat Lactoylglutathione Lyase, GLO1 ELISA Kit
Reactivity	Rat
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	GLO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,820
Uniprot ID	Q6P7Q4
Solubility (25°C)	240ng/ml
Sensitivity	0.29ng/ml
Entrez	294320
Expiration Date	Please enquire.