

Product Datasheet

Rat Steroid 5-alpha-reductase 2 ELISA Kit (orb1211591)

Description	Rat Steroid 5-alpha-reductase 2, SRD5A2 ELISA Kit
Reactivity	Rat
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	SRD5A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,780
Uniprot ID	P24008
Solubility (25°C)	32ng/ml
Sensitivity	0.024ng/ml
Entrez	24950
Expiration Date	Please enquire.