

Product Datasheet

Rat Glyceraldehyde-3-phosphate Dehydrogenase ELISA Kit (orb1211603)

Description	Rat Glyceraldehyde-3-phosphate Dehydrogenase, GAPDH ELISA Kit
Reactivity	Rat
Range	15-3800ng/L
Tested Applications	ELISA
Target	GAPDH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,828
Uniprot ID	P04797
Solubility (25°C)	4000ng/L
Sensitivity	7.42ng/L
Entrez	10835113724383
Expiration Date	Please enquire.