

Product Datasheet

Rat Caspase-2 ELISA Kit (orb1211609)

Description	Rat Caspase-2, CASP2 ELISA Kit
Reactivity	Rat
Range	1-400ng/L
Tested Applications	ELISA
Target	CASP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,728
Uniprot ID	P55215
Solubility (25°C)	480ng/L
Sensitivity	0.52ng/L
Entrez	64314
Expiration Date	Please enquire.