

Product Datasheet

Rat Low Density Lipoprotein Receptor ELISA Kit (orb1211635)

Description	Rat Low Density Lipoprotein Receptor, LDLR ELISA Kit
Reactivity	Rat
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	LDLR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	96,622
Uniprot ID	P35952
Solubility (25°C)	64ng/ml
Sensitivity	0.09ng/ml
Entrez	300438
Expiration Date	Please enquire.