

Product Datasheet

Rat Defensin Beta 2 ELISA Kit (orb1211638)

Description	Rat Defensin Beta 2, DEFB2 ELISA Kit
Reactivity	Rat
Range	5-2000ng/L
Tested Applications	ELISA
Target	DEFB4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	6,946
Uniprot ID	O88514
Solubility (25°C)	2400ng/L
Sensitivity	2.61ng/L
Entrez	64389
Expiration Date	Please enquire.