

Product Datasheet

Rat Calpain-1 Catalytic Subunit ELISA Kit (orb1211640)

Description	Rat Calpain-1 Catalytic Subunit, CAPN1 ELISA Kit
Reactivity	Rat
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	CAPN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	82,119
Uniprot ID	P97571
Solubility (25°C)	40ng/ml
Sensitivity	0.028ng/ml
Entrez	29153
Expiration Date	Please enquire.