

Product Datasheet

Rat Progranulin ELISA Kit (orb1211661)

Description	Rat Progranulin, PGRN ELISA Kit
Reactivity	Rat
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	GRN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,370
Uniprot ID	P23785
Solubility (25°C)	240ng/ml
Sensitivity	0.26ng/ml
Entrez	29143
Expiration Date	Please enquire.