

Product Datasheet

Rat Cxc-Chemokine Receptor 4 ELISA Kit (orb1211696)

Description	Rat Cxc-Chemokine Receptor 4, CXCR4 ELISA Kit
Reactivity	Rat
Range	0.1-30ng/mL
Tested Applications	ELISA
Target	CXCL4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,287
Uniprot ID	P06765
Solubility (25°C)	32ng/mL
Sensitivity	0.06ng/mL
Entrez	360918
Expiration Date	Please enquire.