

Product Datasheet

Rat Acetyl-CoA Carboxylase 1 ELISA Kit (orb1211707)

Description	Rat Acetyl-CoA Carboxylase 1, ACC1 ELISA Kit
Reactivity	Rat
Range	0.05-40ng/ml
Tested Applications	ELISA
Target	ACACA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	265,194
Uniprot ID	P11497
Solubility (25°C)	48ng/ml
Sensitivity	0.03ng/ml
Entrez	60581
Expiration Date	Please enquire.