

Product Datasheet

Rat Endothelial Cell-specific Molecule 1 ELISA Kit (orb1211709)

Description	Rat Endothelial Cell-specific Molecule 1, ESM-1 ELISA Kit
Reactivity	Rat
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	ESM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,075
Uniprot ID	P97682
Solubility (25°C)	96ng/mL
Sensitivity	0.163ng/mL
Entrez	64536
Expiration Date	Please enquire.