

Product Datasheet

Rat Hormone Sensitive Lipase ELISA Kit (orb1211716)

Description	Rat Hormone Sensitive Lipase, HSL ELISA Kit
Reactivity	Rat
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	LIPE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	116,812
Uniprot ID	P15304
Solubility (25°C)	96ng/ml
Sensitivity	0.15ng/ml
Entrez	25330
Expiration Date	Please enquire.