

Product Datasheet

Rat Glucose-6-phosphate 1-dehydrogenase ELISA Kit (orb1211721)

Description	Rat Glucose-6-phosphate 1-dehydrogenase, G6PD ELISA Kit
Reactivity	Rat
Range	0.1-40mU/mL
Tested Applications	ELISA
Target	G6PDX
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,376
Uniprot ID	P05370
Solubility (25°C)	48mU/mL
Sensitivity	0.054mU/mL
Entrez	24377
Expiration Date	Please enquire.