

Product Datasheet

Rat Glypican-3 ELISA Kit (orb1211729)

Description	Rat Glypican-3, GPC-3 ELISA Kit
Reactivity	Rat
Range	10-2000ng/L
Tested Applications	ELISA
Target	GPC3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	67,049
Uniprot ID	P13265
Solubility (25°C)	2400ng/L
Sensitivity	5.29ng/L
Entrez	25236
Expiration Date	Please enquire.