

Product Datasheet

Rat Catalase ELISA Kit (orb1211736)

Description	Rat Catalase, CAT ELISA Kit
Reactivity	Rat
Range	1-300ng/ml
Tested Applications	ELISA
Target	CAT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,757
Uniprot ID	P04762
Solubility (25°C)	320ng/ml
Sensitivity	0.52ng/ml
Entrez	24248
Expiration Date	Please enquire.