

Product Datasheet

Rat Tumor Necrosis Factor A ELISA Kit (orb1211747)

Description	Rat Tumor Necrosis Factor A, TNF-A ELISA Kit
Reactivity	Rat
Range	5-1000ng/L
Tested Applications	ELISA
Target	TNFA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,806
Uniprot ID	P16599
Solubility (25°C)	1280ng/L
Sensitivity	2.51ng/L
Entrez	10369438024835
Expiration Date	Please enquire.