

Product Datasheet

Rat Acetylcholinesterase ELISA Kit (orb1211753)

Description	Rat Acetylcholinesterase, ACHE ELISA Kit
Reactivity	Rat
Range	0.5-150U/L
Tested Applications	ELISA
Target	ACHE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,196
Uniprot ID	P37136
Solubility (25°C)	160U/L
Sensitivity	0.25U/L
Entrez	83817
Expiration Date	Please enquire.