

## Product Datasheet

### Rat Heme Oxygenase 1 ELISA Kit (orb1211762)

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Description</b>         | Rat Heme Oxygenase 1,<br>HO-1 ELISA Kit |
| <b>Reactivity</b>          | Rat                                     |
| <b>Range</b>               | 0.05-20ng/ml                            |
| <b>Tested Applications</b> | ELISA                                   |
| <b>Target</b>              | HMOX1                                   |
| <b>Storage</b>             | 2-8°C                                   |
| <b>Note</b>                | For research use only                   |
| <b>Application notes</b>   | Sandwich                                |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates  |
| <b>MW</b>                  | 33,006                                  |
| <b>Uniprot ID</b>          | <b>P06762</b>                           |
| <b>Solubility (25°C)</b>   | 24ng/ml                                 |
| <b>Sensitivity</b>         | 0.023ng/ml                              |
| <b>Entrez</b>              | <b>24451</b>                            |
| <b>Expiration Date</b>     | Please enquire.                         |