

Product Datasheet

Rat Vascular Endothelial Cell Growth Factor D ELISA Kit (orb1211766)

Description	Rat Vascular Endothelial Cell Growth Factor D, VEGF-D ELISA Kit
Reactivity	Rat
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	VEGFD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,113
Uniprot ID	035251
Solubility (25°C)	48ng/mL
Sensitivity	0.053ng/mL
Entrez	360457
Expiration Date	Please enquire.