

Product Datasheet

Rat Cytochrome P450 2E1 ELISA Kit (orb1211773)

Description	Rat Cytochrome P450 2E1, CYP2E1 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	CYP2E1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	56,627
Uniprot ID	P05182
Solubility (25°C)	24ng/mL
Sensitivity	0.019ng/mL
Entrez	25086
Expiration Date	Please enquire.