

## Product Datasheet

### Rat Aquaporin 2 ELISA Kit (orb1211783)

|                            |                                        |
|----------------------------|----------------------------------------|
| <b>Description</b>         | Rat Aquaporin 2, AQP-2 ELISA Kit       |
| <b>Reactivity</b>          | Rat                                    |
| <b>Range</b>               | 3-900ng/mL                             |
| <b>Tested Applications</b> | ELISA                                  |
| <b>Target</b>              | AQP2                                   |
| <b>Storage</b>             | 2-8°C                                  |
| <b>Note</b>                | For research use only                  |
| <b>Application notes</b>   | Sandwich                               |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates |
| <b>MW</b>                  | 28,931                                 |
| <b>Uniprot ID</b>          | <b>P34080</b>                          |
| <b>Solubility (25°C)</b>   | 960ng/mL                               |
| <b>Sensitivity</b>         | 1.56ng/mL                              |
| <b>Entrez</b>              | <b>25386</b>                           |
| <b>Expiration Date</b>     | Please enquire.                        |