

Product Datasheet

Rat Heat Shock Protein 60 ELISA Kit (orb1211794)

Description	Rat Heat Shock Protein 60, HSP60 ELISA Kit
Reactivity	Rat
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	HSPD1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	60,955
Uniprot ID	P63039
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	63868
Expiration Date	Please enquire.