

Product Datasheet

Rat Hyaluronan Synthase 2 ELISA Kit (orb1211795)

Description	Rat Hyaluronan Synthase 2, HAS2 ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	HAS2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,534
Uniprot ID	035776
Solubility (25°C)	1600ng/L
Sensitivity	8.01ng/L
Entrez	25694
Expiration Date	Please enquire.