

Product Datasheet

Rat Coagulation Factor 10 ELISA Kit (orb1211800)

Description	Rat Coagulation Factor 10, F10 ELISA Kit
Reactivity	Rat
Range	0.5-200ng/mL
Tested Applications	ELISA
Target	F10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,265
Uniprot ID	Q63207
Solubility (25°C)	240ng/mL
Sensitivity	0.27ng/mL
Entrez	29243
Expiration Date	Please enquire.