

Product Datasheet

Rat Secretin ELISA Kit (orb1211847)

Description	Rat Secretin, SCT ELISA Kit
Reactivity	Rat
Range	1-400ng/L
Tested Applications	ELISA
Target	SCT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,072
Uniprot ID	P11384
Solubility (25°C)	480ng/L
Sensitivity	0.52ng/L
Entrez	24769
Expiration Date	Please enquire.