

Product Datasheet

Rat Interleukin 15 ELISA Kit (orb1211860)

Description	Rat Interleukin 15, IL-15 ELISA Kit
Reactivity	Rat
Range	0.5-200ng/L
Tested Applications	ELISA
Target	IL15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,804
Uniprot ID	P97604
Solubility (25°C)	240ng/L
Sensitivity	0.32ng/L
Entrez	25670
Expiration Date	Please enquire.