

Product Datasheet

Rat Mothers Against Decapentaplegic Homolog 7 ELISA Kit (orb1211868)

Description	Rat Mothers Against Decapentaplegic Homolog 7, SMAD7 ELISA Kit
Reactivity	Rat
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	SMAD7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46458
Uniprot ID	O88406
Solubility (25°C)	16ng/ml
Sensitivity	0.025ng/ml
Entrez	81516
Expiration Date	Please enquire.