

Product Datasheet

Rabbit Apolipoprotein A1 ELISA Kit (orb1211883)

Description	Rabbit Apolipoprotein A1, APO-A1 ELISA Kit
Reactivity	Rabbit
Range	2-700ug/ml
Tested Applications	ELISA
Target	APOA1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,591
Uniprot ID	P09809
Solubility (25°C)	800µg/ml
Sensitivity	1.21µg/ml
Entrez	100009253
Expiration Date	Please enquire.