

Product Datasheet

Rabbit Insulin-like Growth Factors 1 ELISA Kit (orb1211884)

Description	Rabbit Insulin-like Growth Factors 1, IGF-1 ELISA Kit
Reactivity	Rabbit
Range	5-1000ng/ml
Tested Applications	ELISA
Target	IGF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,092
Uniprot ID	Q95222
Solubility (25°C)	1280ng/ml
Sensitivity	2.25ng/ml
Entrez	100008668
Expiration Date	Please enquire.