

Product Datasheet

Rabbit Endothelin 1 ELISA Kit (orb1211891)

Description	Rabbit Endothelin 1, ET-1 ELISA Kit
Reactivity	Rabbit
Range	3-900ng/L
Tested Applications	ELISA
Target	EDN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,828
Uniprot ID	P29560
Solubility (25°C)	960ng/L
Sensitivity	1.63ng/L
Entrez	100009270
Expiration Date	Please enquire.