

Product Datasheet

Rabbit Tumor Necrosis Factor A ELISA Kit (orb1211893)

Description	Rabbit Tumor Necrosis Factor A, TNF-A ELISA Kit
Reactivity	Rabbit
Range	0.5-200ng/L
Tested Applications	ELISA
Target	TNFA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25817
Uniprot ID	P04924
Solubility (25°C)	240ng/L
Sensitivity	0.24ng/L
Entrez	100009088
Expiration Date	Please enquire.