

Product Datasheet

Rabbit Angiotensin Converting Enzyme ELISA Kit (orb1211897)

Description	Rabbit Angiotensin Converting Enzyme, ACE ELISA Kit
Reactivity	Rabbit
Range	0.5-180ng/ml
Tested Applications	ELISA
Target	ACE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	150,406
Uniprot ID	P12822
Solubility (25°C)	200ng/ml
Sensitivity	0.33ng/ml
Entrez	100009274
Expiration Date	Please enquire.