

Product Datasheet

Rabbit C-Reactive Protein ELISA Kit (orb1211898)

Description	Rabbit C-Reactive Protein, CRP ELISA Kit
Reactivity	Rabbit
Range	20-4500ng/L
Tested Applications	ELISA
Target	CRP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,493
Uniprot ID	P02742
Solubility (25°C)	4800ng/L
Sensitivity	9.91ng/L
Entrez	100009091
Expiration Date	Please enquire.