

Product Datasheet

Porcine Vascular Endothelial Growth Factor A ELISA Kit (orb1211933)

Description	Porcine Vascular Endothelial Growth Factor A, VEGFA ELISA Kit
Reactivity	Porcine
Range	10-3000ng/L
Tested Applications	ELISA
Target	VEGFA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,368
Uniprot ID	P49151
Solubility (25°C)	3200ng/L
Sensitivity	5.19ng/L
Entrez	397157
Expiration Date	Please enquire.