

Product Datasheet

Porcine Ghrelin ELISA Kit (orb1211943)

Description	Porcine Ghrelin, GHRL ELISA Kit
Reactivity	Porcine
Range	3-700ng/L
Tested Applications	ELISA
Target	GHRL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,786
Uniprot ID	Q9GKY5
Solubility (25°C)	800ng/L
Sensitivity	1.52ng/L
Entrez	396728
Expiration Date	Please enquire.