

Product Datasheet

Porcine Transthyretin ELISA Kit (orb1211945)

Description	Porcine Transthyretin, TTR ELISA Kit
Reactivity	Porcine
Range	5-2000ng/L
Tested Applications	ELISA
Target	TTR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,081
Uniprot ID	P50390
Solubility (25°C)	2400ng/L
Sensitivity	2.53ng/L
Entrez	397419
Expiration Date	Please enquire.