

Product Datasheet

Porcine Vascular Cell Adhesion Molecule 1 ELISA Kit (orb1211955)

Description	Porcine Vascular Cell Adhesion Molecule 1, VCAM-1 ELISA Kit
Reactivity	Porcine
Range	10-2000ng/mL
Tested Applications	ELISA
Target	VEGFA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,368
Uniprot ID	P49151
Solubility (25°C)	2400ng/mL
Sensitivity	5.35ng/mL
Entrez	397157
Expiration Date	Please enquire.