

Product Datasheet

Porcine Endothelin 1 ELISA Kit (orb1211958)

Description	Porcine Endothelin 1, ET-1 ELISA Kit
Reactivity	Porcine
Range	3-700ng/ml
Tested Applications	ELISA
Target	EDN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,259
Uniprot ID	P09558
Solubility (25°C)	800ng/ml
Sensitivity	1.52ng/ml
Entrez	396915
Expiration Date	Please enquire.